EV-Perm™ Permeabilization Pretreatment Kit for Exosome Membrane

【Product Information】

Code No.	Product Name	Volume	Storage
294-85701	EV-Perm™ Permeabilization Pretreatment Kit for Exosome Membrane	1 kit	Frozen (-20℃)

【Kit Contents】

Reagent	Volume	Storage
Reagent A (20×)	500 μL×1	Frozen (-20℃)
Reagent B (10×)	5 mL×1	Before thawing : -20℃ After thawing : 2~10℃
Reagent C (100×)	500 μL×1	Before thawing : -20℃ After thawing : 2~10℃

【How to store each reagent】

Store each reagent according to the storage condition above.

Preparation of reagents

1. Reagent A (20×)

Following thaw, mix the vial by gently tapping, and prepare aliquots. Close the lid tightly and immediately store the remaining solution in the freezer (-20℃).

- If stored in a liquid state for a long time, components may precipitate.
- If precipitation is observed, heat it at 60℃ for about 30 minutes to dissolve it, then store it in the freezer.

2. Reagent B (10×) and Reagent C (100×)

Following thaw, mix the vial by gently tapping, and prepare aliquots. Close the lid tightly and immediately store the remaining solution in the refrigerator (2~10℃).

[About this product]

This product is a reagent that enables the detection of internal markers of extracellular vesicles. Until now, the detection targets of the PS Capture™ Exosome ELISA Kit Series and PS Capture™ Exosome Flow Cytometry Kit were limited to the detection of surface markers of extracellular vesicles. However, when used in combination with this product, it is possible to improve the permeability of the extracellular vesicle membrane surface and detect internal markers.

For PS Capture™ Exosome ELISA Kit

[Equipment and reagents to be supplied by user]

① ELISA Kit (Check) :

- One of the following kits

Code No.	Product Name
297-79201	PS Capture™ Exosome ELISA Kit (Anti Mouse IgG POD)
298-80601	PS Capture™ Exosome ELISA Kit (Streptavidin HRP)

② Reagents (Check) :

- Distilled Water
 Specific antibody for Extracellular vesicle internal protein

If using the PS Capture™ Exosome ELISA Kit (anti-mouse IgG POD), please prepare a mouse monoclonal antibody, and if using the PS Capture™ Exosome ELISA Kit (Streptavidin HRP), prepare a biotin-labeled antibody.

③ PS Capture™ Exosome ELISA Kit included reagents

- Washing Buffer (10×)
 Exosome Binding Enhancer (100×)
 Reaction Buffer

Please use the Exosome Binding Enhancer and Reaction Buffer in the PS Capture™ Exosome ELISA Kit. Even if you use Buffer with this kit, you will still have Buffer for the PS Capture™ Exosome ELISA Kit.

④ Equipments (Check) :

- Microwell plate (PS Capture™ Exosome ELISA kit accessory)
 Plate seal (PS Capture™ Exosome ELISA kit accessory)
 Microplate shaker
 Microplate reader

【Preparation of samples】

For PS Capture™ Exosome ELISA Kit (anti-mouse IgG POD)

It is necessary to optimize the concentration of extracellular vesicles used to detect the target marker using this product. The optimal concentration of extracellular vesicles was determined using the PS Capture™ Exosome ELISA kit (anti-mouse IgG POD) such that the resulting absorbance (450 nm minus 620 nm) was in the range of 0.2 to 2.5. Prepare. To adjust the concentration of extracellular vesicles, dilute them with the reaction/wash solution (1x) included in the kit.

For PS Capture™ Exosome ELISA Kit (Streptavidin HRP)

It is necessary to optimize the concentration of extracellular vesicles used to detect the target marker using this product. The optimal concentration of extracellular vesicles was measured using the PS Capture™ Exosome ELISA Kit (Streptavidin HRP) and adjusted so that the resulting absorbance (450 nm minus 620 nm) was in the range of 0.2 to 2.5. To do. To prepare the concentration of extracellular vesicles, dilute at least two times with the Reaction Buffer included in the kit.

【Preparation of reagents】

1. Preparation of permeabilization solution

Dilute each reagent, referring to Table 1 below. After dilution, mix well. Table 1. Dilution ratio and amount of reagent added for each reagent used in permeabilization solution

Reagent	Dilution	12 samples* (24 wells)	24 samples* (48 wells)	48 samples* (96 wells)
Distilled Water	-	2.1 mL	4.2 mL	8.4 mL
Reagent A (20×)	1:20	125 µL	250 µL	500 µL
Reagent B (10×)	1:10	250 µL	500 µL	1.0 mL
Reagent C (100×)	1:100	25 µL	50 µL	100 µL
Total	-	2.5 mL	5.0 mL	10.0 mL

※ For performing in duplicate.

2. Preparation of negative control solution

Dilute each reagent referring to Table 2 below. After diluting, mix well

(Do not add Reagent A when preparing the negative control solution.)

Table 2. Dilution ratio and amount of reagent added for each reagent used for negative control

Reagent	Dilution	12 samples* (24 wells)	24 samples* (48 wells)	48 samples* (96 wells)
Distilled Water	-	2.2 mL	4.5 mL	8.9 mL
Reagent B (10×)	1:10	250 µL	500 µL	1.0 mL
Reagent C (100×)	1:100	25 µL	50 µL	100 µL
Total	-	2.5 mL	5.0 mL	10.0 mL

※ For performing in duplicate.

3. Preparation of washing solution (1x)

Washing Buffer (10x) included with the PS Capture™ Exosome ELISA Kit was diluted 10 times with purified water. After that, add 1/100 amount of Exosome Binding Enhancer (100x) to the diluted solution and mix well. 10

【Procedure】

1. Process Flow

Process		Reagent · Equipment	Time
1.	Capture Exosome	Washing the plate	Washing buffer (1×)
		Capture Exosome on the plate	Sample Capture (1× Reaction Buffer)
		Washing the plate	Washing buffer (1×)
2.	Permeabilization	Permeabilization*	① permeabilization solution ② negative control solution
		Washing the plate	Washing buffer (1×)

※ The permeabilization solution and negative control solution are separate for each sample and detection condition in the permeabilization process.

The following steps are performed according to the PS Capture™ Exosome ELISA kit.

Process		Reagent · Equipment	Time
3.	1st antibody reaction	1st antibody reaction	Add detection antibody
		Washing the plate	Washing buffer (1×)
4.	2nd antibody / streptavidin-HRP reaction	2nd antibody / streptavidin-HRP reaction	2nd antibody / streptavidin-HRP reaction
		Washing the plate	Washing buffer (1×)
5.	Substrate	Substrate addition	TMB Solution Stop Solution
		ELISA plate reading	microplate reader

2. The procedure of #2 permeabilization in the process flow above.

- (1) Follow the instructions for the PS Capture™ Exosome ELISA kit and proceed to the washing step Procedure# 4 before adding the 1st antibody reaction solution.
- (2) Add 100 µL of permeabilization solution or 100 µL of negative control solution to the extracellular vesicle sample (well).

- (3) Apply a plate seal and incubate at room temperature for 1 hour while stirring at 500 rpm using a microplate shaker.
- (4) Decant the reaction solution and wash each well three times with 300-350 μ L of washing solution (1x).
- (5) From the process flow "3. 1st antibody reaction", proceed to step "[Procedure] 5)" according to the PS Capture™ Exosome ELISA Kit instruction manual.

Please refer to the PS Capture™ Exosome ELISA Kit Instruction Manual for post-measurement calculations and analysis methods.

For PS Capture™ Exosome Flow Cytometry Kit

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