

## SPECIFICATION

## Sodium Bisulfite

Japanese Pharmacopoeia

REQUIREMENT	SPECIFICATION
* Description	White shot or powder
Identification	to pass test
Purity	-
(1) Clarity and color of solution	to pass test
(2) Thiosulfate	to pass test
(3) Iron	not more than 20ppm
(4) Arsenic	not more than 4ppm
Assay (as SO <sub>2</sub> )	64.0~67.4%

Japanese Pharmacopoeia(\*:Additional test performed by Wako)